



## SOP Title: qPCR for *Brucella* spp.

SOP ID: 309, Version: 1.0

Effective Date: 3 March, 2026

Review Date: 3 March, 2027

Primary Author (s): Tom Edwards, Caitlin Greenland-Bews

Secondary Author: Eric Fevre

Approver: Siobhan Mor

### 1. Purpose

To ensure a standardized and consistent approach to the PCR-based detection of *Brucella* spp. DNA in whole blood samples collected as part of the ZAFI project.

### 2. Scope

This SOP applies to all laboratory personnel involved in the processing, analysis, and interpretation of samples for *Brucella* spp. PCR testing under the ZAFI project across Kenya, Uganda, and Ethiopia.

### 3. Abbreviations

- **qPCR:** Quantitative Polymerase Chain Reaction

### 4. Responsibilities

5. Role	Responsibility
Laboratory technician	Perform laboratory testing according to this SOP; handle and process samples appropriately; adhere to biosafety and quality control procedures; document results accurately in logbooks and KOBO tools; report deviations or issues to the laboratory lead/supervisor.
Laboratory lead/supervisor	Oversee laboratory testing; ensure staff competency and adherence to quality and biosafety standards; review and validate results; escalate issues as appropriate.
Site PI	Ensure relevant staff are trained on this SOP.

### 6. Safety and Precautions

- Preparation of mastermix should be conducted in a “Clean room” or “Clean hood” that is separate to template DNA and any amplification procedures to prevent contamination
- Template (DNA) and samples should be stored separately to mastermix and assay reagents
- When calculating the volume of mastermix to make, include enough volume for a single positive control per target and a negative control

## 7. Consumables and Reagents

Item	
<b>Reagents</b>	<b>Equipment</b>
TaqMan Fast Advanced Master Mix (2x)	Microcentrifuge
Primers and Probes	Pipettes (10µl, 20 µl, 200 µl, 1000 µl)
Nuclease Free Water	Vortex
Bcsp31_std positive control material (gBlock)	<b>Tubes/Plates/other</b>
IS711_std positive control material (gBlock)	Plastics consumables
	Pipette tips
	Disposable gloves

## 8. qPCR Procedure

### 7.1 Preparing the PCR Reaction Mix

- Retrieve the working (10uM) stocks of primers, probes and Master mix for the *Brucella* qPCR (Table 1) **including the 0.3uM stocks of RNaseP primer and probe.**
- When the reagents are completely thawed vortex each tube briefly and spin in microcentrifuge to collect liquid at the bottom.
- Calculate the number of reactions needed for the samples (one reaction per sample) **including a positive control for each target and a negative control.** Multiply the volume per reaction of each component by the total number of reactions.
- Combine these volumes into one 1.5ml Eppendorf **\*This should be done in a separate room/hood to template addition.\***

Table of reactions			
Component	Final concentration	Volume per Reaction	Total volume to add to the 1.5ml Eppendorf assuming <i>n</i> reactions
		96-well (0.1ml) plates	
<b>Reaction mix</b>			
TaqMan Fast Advanced Master Mix (2x)	1X	5.50	5.50 x <i>n</i>
bcsp31_1163_F (10uM)	300nM	0.42	0.42 x <i>n</i>
bcsp31_1163_R (10uM)	300nM	0.42	0.42 x <i>n</i>
bcsp31_Tq_probe-FAM (10uM)	200nM	0.28	0.28 x <i>n</i>
IS711_421_F (10uM)	300nM	0.42	0.42 x <i>n</i>
IS711_Tq_R_ (10uM)	300nM	0.42	0.42 x <i>n</i>
IS711 probe-VIC (10uM)	50nM	0.28	0.28 x <i>n</i>
RNase P primer mix (0.3uM)		0.42	0.42 x <i>n</i>
RNase P probe mix (TEXASRED (0.3uM))		0.17	0.17 x <i>n</i>
NFW	-	2.7	2.7 x <i>n</i>
<b>Aliquot per well</b>			
Reaction mix		10ul	
DNA		2.5ul	
<b>Total volume</b>		<b>12.5ul</b>	

- Vortex briefly to mix and spin to ensure all liquid is at the bottom.
- Centrifuge in a microcentrifuge to bring reaction mix to the bottom of the tube.
- Aliquot 10µl volume of PCR Reaction Mix to each well of the plate

- h) Retrieve DNA eluates immediately prior to template addition and thaw on ice
- i) Transfer plate to template addition room/space and add 2.5µl DNA, or nuclease-free water for Non-Template Control, or positive control to separate wells.
- j) Seal the reaction plate with optical adhesive film or strip-caps and then centrifuge briefly to bring the PCR reaction mix to the bottom of the well.

### 8.2. Setting up the machine

- a) Set up the following run profile on the qPCR machine.
- b) Ensure that each well is set to collect data in the FAM, VIC and TexasRed channels

Step	Temperature	Time	Cycles	Fluorescence acquisition (channel)
Hold	50	2:00	1	
Hold	95	0:20	1	
Denature	95	0:03	40	
Anneal/Extend	60	0:30		FAM, VIC and TEXASRED

### 7.3 Interpretation

- a) The negative control must have no amplification present in any channel
- b) The positive controls must have amplification in their respective channels
  - Bcsp standard (positive control) should amplify in the FAM channel
  - IS711 standard (positive control) should amplify in the VIC Channel
- c) Positive amplification should be observed in the TexasRed channel for all human samples
- d) Results with a Ct below 38 in either the FAM or VIC channel can be considered as positive.

## 9. Record Keeping

- Enter results into the results sheet
- Enter results into the dedicated **KOBO form**

## 10. Primer Sequences

Primer/probe ID	Sequence
bcp31_1163_F	TCTTTGTGGGCGGCTATCC
bcp31_1199_R	CCGTTGAGATGGCCAGTT
bcp31_Tq (probe)	FAM-ACGGGCGCAATCT-MGB-NFQ
IS711_421_F	CGCTCGCGCGGTGGAT
IS_Tq_R	CTTGAAGCTTGCAGACAGTCACC
IS711_probe	VIC-ACGACCAAGCTGCATGCTGTTGTGCGATG-BHQ1
RNaseP_F	CCAAGTGTGAGGGCTGAAAAG
RNaseP_R	TGTTGTGGCTGATGAACTATAAAAAGG
RNaseP_probe	TEXASRED-CCCCAGTCTCTGTGAGCACTCCCTTC-BHQ

## 11. Review and Approval

Approved By: *Prof Siobhan Mor*

Title: *Chief Investigator, ZAFI*

Date: *28 February 2026*