

## SOP: IgG IFA for *Coxiella burnetii* (Q fever)

SOP ID: 322, Version: 1.1

Effective Date: 15 April, 2026

Review Date: 15 April, 2027

Primary Author (s): Tom Edwards, Caitlin Greenland-Bews

Secondary Author: Eric Fevre

Approver: Siobhan Mor

### 1. Purpose

To describe the standardized procedure for detecting **IgG antibodies against *Coxiella burnetii* in human serum using the Indirect Immunofluorescent Assay (IFA)**, following the manufacturer's instructions and as adapted to the ZAFI project.

### 2. Scope

This SOP applies to all laboratory personnel performing and interpreting the results of the *Coxiella burnetii* IFA test at ZAFI project laboratories in Uganda, Kenya, and Ethiopia.

### 4. Abbreviations

- **IFA:** Indirect Immunofluorescent Assay
- **IgG:** Immunoglobulin G
- **PBS:** Phosphate Buffered Saline
- **QC:** Quality Control
- **FITC:** Fluorescein Isothiocyanate

### 3. Responsibilities

Role	Responsibility
<b>Laboratory technician</b>	Perform laboratory testing according to this SOP; handle and process samples appropriately; adhere to biosafety and quality control procedures; document results accurately in logbooks and KOBO tools; report deviations or issues to the laboratory lead/supervisor.
<b>Laboratory lead/supervisor</b>	Oversee laboratory testing; ensure staff competency and adherence to quality and biosafety standards; review and validate results; escalate issues as appropriate.
<b>Site PI</b>	Ensure relevant staff are trained on this SOP.

## 5. Materials and Equipment

Item	
<b>Supplied in the Kit:</b>	<b>Additional Materials (not supplied):</b>
<i>Coxiella burnetii</i> IFA slides (10 slides × 10 wells, Phase II antigen)	Calibrated micropipettes and tips
PBS powder (pH 7.2)	Humid incubation chamber
Positive control serum (IgG)	Thermostatically controlled incubator (37 °C)
Negative control serum	Fluorescence microscope (400× magnification) with appropriate filters
FITC-labelled anti-human IgG conjugate	Distilled or deionized water
Mounting medium (buffered glycerol)	Coverslips (24 × 60 mm)
	Timer
	PPE: lab coat, gloves, eye protection

## 6. Safety Precautions

- Handle all human specimens as potentially infectious
- Wear appropriate PPE at all times
- Avoid skin and eye contact with reagents containing sodium azide and Evans blue
- Dispose of waste according to the SOP on Biohazard Waste Disposal
- Slides contain inactivated antigen but must still be handled with care

## 7. Procedure

### 7.1 Preparation

- a) Bring all samples and reagents to room temperature before use and allow the slides to reach room temperature before opening.
- b) Reconstitute PBS by dissolving the contents of one vial in 1 litre of distilled water. Mix thoroughly until fully dissolved
  - Store reconstituted PBS at **2–8 °C** (maximum 4 months; never beyond expiry date)

### 7.2 Sample Dilution

Initial screening is performed on 1/64 and 1/128 dilutions:

- a) Prepare a **1/64 dilution** of the serum by adding 10 µL sample + 630 µL PBS
- b) Prepare **1/128 dilution** by two-fold dilution of the 1/64 dilution using 50ul PBS and 50ul of the 1/64 diluted sample

**NOTE: DO NOT DILUTE CONTROL SERA (VIRCELL COXIELLA IgG POSITIVE CONTROL and VIRCELL COXIELLA NEGATIVE CONTROL)**

### 7.3 Slide Preparation and Incubation

- a) Apply **20 µL** of diluted samples (1/64 and 1/128) in two slide wells
- b) Add **20 µL** of positive and negative controls to separate wells
- c) Incubate slides in a humid chamber at **37 °C for 30 minutes**

### 7.4 Washing

- a) Rinse slides gently with PBS (avoid directing stream into wells)
- b) Immerse slides in PBS for **10 minutes**
- c) Briefly dip slides in distilled water
- d) Allow slides to air-dry completely

### 7.5 Conjugate Incubation

- a) Add **20 µL** of FITC-labelled anti-human IgG conjugate to each well (No dilution required)
- b) Incubate at **37 °C for 30 minutes** in a humid chamber
- c) Rinse slides gently with PBS (avoid directing stream into wells)
- d) Immerse slides in PBS for **10 minutes**
- e) Briefly dip slides in distilled water
- f) Allow slides to air-dry completely

### 7.6 Mounting and Reading

- a) Add a small drop of mounting medium to each well
- b) Carefully apply coverslip
- c) Read slides **immediately** using a fluorescence microscope at **400× magnification**
- d) If not read immediately, store slides at **2–8 °C in the dark** for up to **24 hours**

### 7.7 Further Dilutions of Positive Screens

If the screening dilutions are positive, further analyse with increasing dilutions up to 1/2048:

- a) Prepare two-fold dilutions up to 1/2048 starting from the already prepared 1/128 dilution.
  - i. Add 50ul PBS and 50ul of the 1/128 diluted sample to prepare a 1/256 dilution.
  - ii. Add 50ul PBS and 50ul of the 1/256 dilution to prepare a 1/512 dilution
  - iii. Add 50ul PBS and 50ul of the 1/512 dilution to prepare a 1/1024 dilution
  - iv. Add 50ul PBS and 50ul of the 1/1024 dilution to prepare a 1/2048 dilution
- b) Analyse according to Section 7.6

## 8. Interpretation of Results

Result	Interpretation
Apple-green fluorescence with cocco-bacillary morphology	Positive
No specific fluorescence	Negative

- The **antibody titre** is the highest dilution showing positive fluorescence
- Titres **≥1/64** are generally considered indicative of infection

- Seroconversion or a significant rise in titre between paired samples confirms recent infection
- The reaction is negative when no fluorescence can be observed

### 9. Quality Control

- Include **positive and negative controls** in every run
- Results are valid only if:
  - Positive control shows characteristic fluorescence
  - Negative control shows no fluorescence
- Failed QC invalidates the run and requires repeat testing

### 10. Documentation

- Record assay details in the **IFA Run Log**, including date, operator, kit lot number, sample IDs, dilutions, final titres and interpretation
- Enter validated results into dedicated **KOBO form**
- Maintain **IFA QC logs** for review.
- Archive raw data printouts and electronic files securely.

### 11. Related templates

- IFA Run Log
- IFA QC Log

### 12. Review and Approval

Approved By: *Prof Siobhan Mor*

Title: *Chief Investigator, ZAFI*

Date: *15 April 2026*

## Annex

### IFA Slide Layout Template




Figure 1: Diagram showing an example slide layout. PC = Positive Control, NC= Negative Control. The two dilutions are plated next to each other. S1 and S2 are names used to demonstrate two different samples, please use the sample numbers and keep a note of the well positions of each sample.